

Pherocious™ CT

Apheresis Catheter

apheresis with power and vision

A staggered medial/proximal lumen location and power injection capability at the distal tip make the Pherocious™ CT catheter the new standard for addressing the challenges of long-term apheresis.

The first power injectable apheresis catheter, for patient comfort and needle-free contrast media injection

10.5F and 12.5F sizes with high capacity lumens for the broad range of apheresis uses

- 10.5F - up to 80ml/min flow rate
- 12.5F - up to 150ml/min flow rate

Atraumatic tip

Power injection
distal lumen

Staggered medial and
proximal lumens



r4 Vascular, Inc.
7550 Meridian Circle N, Suite 150
Maple Grove, MN 55369
1-866-943-8090

www.r4vascular.com

Pherocious™ CT

Apheresis Catheter

Staggered aspiration and infusion positions, with power injection capability - all in one triple lumen catheter



The Pherocious™ CT Apheresis Catheter brings power injection with distinct infusion and aspiration zones to long-term access, designed for increased patient comfort and more reliable long-term use.

Order#	Description	French Size	Length		Lumen Gauge power/medial/proximal
			hub to tip (cm)	cuff to tip (cm)	
43102440	Apheresis Catheter Kit	10.5	24	19	17/ 16/ 15
43102840	Apheresis Catheter Kit	10.5	28	23	17/ 16/ 15
43103240	Apheresis Catheter Kit	10.5	32	27	17/ 16/ 15
43122440	Apheresis Catheter Kit	12.5	24	19	17/ 15/ 12
43122840	Apheresis Catheter Kit	12.5	28	23	17/ 15/ 12
43123240	Apheresis Catheter Kit	12.5	32	27	17/ 15/ 12

All sizes rated for 5ml/sec at 300psi maximum power injection rate. One catheter per pack.

Kits contain:

- Pherocious CT catheter
- valved tear-away introducer (12.5F Kits only)
- tunneler
- 0.035" X 70cm J-tip guidewire
- three apheresis Luer caps
- 10ml disposable Luer lock syringe
- dressing
- introducer needle
- 15cm dilator (10.5F sizes - 8F dilator/ 12.5F sizes - 10F dilator)



r4 Vascular, Inc.
7550 Meridian Circle N , Suite 150
Maple Grove, MN 55369
1-866-943-8090
www.r4vascular.com

All products delivered sterile. Single-use devices.
U.S. Federal Law restricts sale to or on the order of a physician.

Pherocious is a trademark of r4 Vascular, Inc.
©2009, r4 Vascular, Inc.